Mufla

Indicação:

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For use in denture polymerization.

Characteristics / Benefits:

- Unprecedented model.
- For use in conventional and microwave techniques.
- Elaborated in polymer of the high rigidity (Engineering Plastic).
- Excellent thermal property.
- Superior mechanical resistance.
- Stainless steel tightening system.
- Ease in deflasking due to self-expelling internal angles.
- It provides ease of inclusion since two parts adapt to each other in a unique way, with no possibility of closing it erroneously (Muffle Base and Counter Muffle).

Presentation:

Unit composed of Muffle Base and Counter Muffle, screws for tightening, multipurpose key and instructions for use.

Composition:

Engineering Plastic and stainless steel inserts.

Instructions for Use:

The MUFLA VIPI STG was developed for the polymerization process of acrylic resins, in the conventional and microwave techniques; Even that the internal walls of MUFLA VIPI STG are ultra-polished, we suggest the application of petroleum jelly in paste for a better deflasking;

Apply vaseline in paste at the borders of the muffle before pressing of the resin, thus avoiding possible excesses impregnated in the muffle.

For the preparation of the Base, proceed in the following manner:

Position the deflasking disc at the opening of the Base;

Pour the plaster on the base of the MUFLA VIPI STG;

Position the mold at the center of the base and complete with plaster up to the limit of the border of the Base.

Observation: Do not leave plaster residue at the border of the Muffle.

After setting of the plaster, insulate it as in the conventional technique.

We recommend, in this phase, the application of the silicon wall [LABOR MASS] for obtainment of better results and greater ease in deflasking.

To pour the plaster complement, proceed in the following manner:

Position the Counter Muffle and tighten the 4 screws of the set;

Prepare a more fluid plaster mass, pouring it through the larger opening of the Counter Muffle.

ACRILIZATION:

Prepare the acrylic resin, in accordance with the manufacturer's orientations, considering the technique used (conventional or microwave), transporting it inwards the Muffle and place a plastic separating sheet on the resin.

Close the MUFLA VIPI STG and take to the press, performing a pressing test.

Open the Muffle, trim the excess, insulate the model, making sure that there is no plaster or resin residue on the borders of the Muffle and openings of the fittings. Close again with the plastic separating sheet and perform the definitive pressing.

After the final pressing (1000 Kgf), and still under pressure install the screws until their washers lie on the muffle's surface, and the Allen key that goes with the product, perform only a quarter (1/4) of a turn after the screw washer leans flat.

Recommendations: For the microwave technique we recommend the use of the oven with rotating plate.

Regardless of the technique used, we recommend the cooling of the Mufila VIPI STG at room temperature.

For opening of the MUFLA VIPI STG just use the appropriate Allen key (flat tip), by inserting it into the lateral pins of the Muffle's body.

For deflasking always use a rubber or plastic hammer.

MUFFLE DIMENSIONS:

External measurement: 12 cm x 12 cm

Internal measurement: 8 cm x 10 cm

Base Height: 24 mm

Counter Muffle Height: 30 mm

IMPORTANT ATTENTION: AFTER THE FINAL TIGHTENING, STILL UNDER PRESSING, PLACE THE SCREWS BY TIGHTENING ONLY A QUARTER (1/4) OF A TURN AFTER IT LIES FLAT WITH THE WASHER.

POLYMERIZATION CYCLE THROUGH MICROWAVE POWER

Indicated for VIPI WAVE STG Resin.

The cycles vary depending on microwave power.

	For 500 W oven		For 800 to 1100W oven		For 1200 to 1400W oven	
INITIAL STAGE	20 min	With 20/30% Power	20 min	With 10/20% Power	20 min	With 10% Power
FINAL STAGE	+ 5 min	With 80/100% Power	+ 5 min	With 50/60% Power	+ 5 min	With 30/40% Power

Muffle Kiln

Warnings / Precautions:

Whenever the MUFLA VIPI STG is used, during pressings, pass a layer of petroleum jelly in paste on the borders of the Muffle, preventing the resin to stick on the manufacturing material during the years, thus giving a longer life to the equipment.

Shelf Life:

Undetermined shelf life, with guarantee of 2 years.

Product declared exempt from registration by the Ministry of Health.

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