ORTO COR ORTO COR MIX

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Resina Acrílica

Indicação:

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Acrylic Resin

Indication:

Self-polymerizable acrylic resin for orthodontic works. ORTO COR is indicated for powder/liquid technique. ORTO COR MIX is indicated for manipulation technique. To replace part of chewing system of the human body.

Contraindications:

There are no known adverse side effects, when the product is used according to usual procedures. Material contains polymerizable monomers which may cause skin sensitization (allergic contact dermatitis) or other allergic reactions in susceptible persons.

Warnings / Precautions:

- Exclusively for dental use;
- Liquid: Flammable, non-toxic material (Flash point: 10°C / 50°F);
- Resin: It can form flammable air/powder mixtures;
- As in all acrylic resins, handling of the product is recommended in ventilated places, using nitrile gloves and protective goggles; however, handling is not recommendable for pregnant women, nursing mothers and people with known allergies to this type of product. It can cause skin irritation in susceptible people;
- Keep away from the light and heat;
- Do not dispose the wastes in sewers and watercourses or in household garbage;
- Product that are not compatible with the resin/liquid: alcohol, water, chloroform, acetone, vaseline and other solvents.
- Dispose the product in accordance with local regulations.
- The powder contains an organic polymerization initiator which may degrade if stored at excessively high temperature.

Instructions for Use:

Prior to starting any work, the containers must be completely clean and free of traces of wax or grease.

In order to obtain better results in the use of ORTO COR resin, it is recommended to SHAKE BEFORE USING.

The basic and essential prerequisites for elaboration of orthodontic apparatuses should be respected. Before starting any work, the parts to be used should be completely clean and free from traces of wax or fat, so that the chemical union between the resins may be perfect. To increase the adherence of the resin on metal parts the use of Silano VIPI-TONE is recommended.

In the indicated applications for the product, good hydration of the model and an adequate insulation is important.

In the powder/liquid technique an amount of powder is placed on the model and saturated with dropped liquid. This way various layers are applied until it reaches the desired thickness. The ORTO COR resin has minimum drainage and for this reason the work can be completed at shorter time.

In the manipulated technique, use ORTO COR MIX, by placing the necessary amount of liquid in the glass jar, silicon or Dappen, in accordance with the amount of acrylic. On the liquid, place the resin until the two are mixed without having excess powder or liquid; Recommended ratio: powder/liquid 2:1 (14g of powder for 7ml of liquid). Homogenize carefully with the spatula, avoiding the incorporation of air bubbles by applying the mass on the model.

After application, eliminate the excess material and, still in the plastic phase, eliminate the bubbles through a hydrostatic press in a polymerizer (place the necessary amount of water, apply 20 pounds of compressed air and allow for polymerization to complete for 20 minutes). After total polymerization, start the finish with bits and mills, then with abrasive stones and conclude with fine textured sandpapers. Perform polishing with soft brush, by using VIPI BRIL ACABAMENTO E POLIMENTO.

ORTO COR resin is especially elaborated, in different colors, to expand its field of use by professionals of the field, aside from motivating its use by patients. Biocompatible pigments have high stability. The resins can be mixed in order to obtain new matrices. High stability of colors and abrasion resistance.

NOTE: Due to the presence of Crosslink, VIPI recommends mechanical polishing only for this product.

Characteristics/Benefits:

- ORTO COR resin is acrylate free, reducing the probability of allergic reactions.
- It supports the most varied polymerization methods.
- Range of colors that assures customization of your work.
- Selected raw materials that assure excellent results in the final product.
- Liquid with crosslink, assuring greater product resistance.
- Rigorous selection and raw material receipt tests combined with successive tests during all manufacturing cycles, under ISO standard requirements, assure the specified results in the final product.

Presentations:

ORTO COR MIX: 80g, 225g, 450g and 2,250kg.

ORTO COR: 80g, 225g, 450g, 2,250kg.

ORTO COR CONCENTRATE: 20g.

LIQUID: 50ml, 120ml, 250ml and 5 liters.

Colors:

NORMAL: Crystal, Blue, Yellow, Green, Ruby and Colorless.

CONCENTRATES: Green, Red, Blue, Black, White and Yellow.

FLUORESCENT: Pink and Lemon Green.

ORTO COR ORTO COR MIX

Technical Datasheet
English

Acrylic Resin

Composition:

ORTO COR MIX
Polymethylmethacrylate
Benzoyl Peroxide
Pigments

ORTO COR
Polymethylmethacrylate
Benzoyl Peroxide
Pigments

LIQUID
Methylmethacrylate
DMT
Inhibitor
EDMA (crosslink)
Fluorescent

Recommendations to patients:

To increase the durability of the prosthesis, daily hygiene is recommended without the use of aggressive substances, such as acids or alkalis.

Shelf Life:

POWDER: 10 years as of the date of manufacture.
liquid: 2 years as of the date of manufacture.

ANVISA Registration nº.:

POWDER: 10216040031
LIQUID: 10216040032

Technical Renposible:

Luiz Paulo Ribeiro Junior - CRQ IV (Regional Council of Chemistry) 03251745

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