### Resina Acrílica Odontológica

**Indicação:**

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### Acrylic Resin

**Indication:**

Colorless and pigmented Self-Polymerizable Acrylic Resin, used in Repairs, Additions, Rebases of Denture and RPD, To replace part of chewing system of the human body.

**Contraindications:**

The material contains polymerizable monomers that can cause skin sensitization (allergic contact dermatitis) or other allergic reactions in susceptible people.

**Warnings/Precautions:**

- For dental use only.
- Liquid: Flammable, non-toxic material (Flash point: 10°C / 50°F).
- Resin: Can form flammable air / dust mixtures.
- Dispose of the product in accordance with local regulations.
- The powder contains an organic polymerization initiator that can degrade if stored at excessively high temperature.
- It is recommended to handle the product in ventilated places, using nitrile gloves and glasses of protection; however, handling is not recommended for pregnant women, nursing mothers and people with recognized allergies to this type of product.
- May cause skin irritation in susceptible people.
- Keep protected from light and heat, closed in a dry and cool place.
- Do not dispose of waste in sewers and watercourses or in domestic waste.
- Products that are not compatible with the resin / monomer: alcohol, water, chloroform, acetone, petroleum jelly and other solvents.
- The base of artificial teeth must be prepared according to the Instructions for Use for the proper adhesion of the prosthesis.
- The use of wax solvent is not recommended. The residual solvent can cause less connection of the teeth to the acrylic base.
- Prohibited to reprocess.

**Step-By-Step Instructions:**

Shake the bottle before using.

**PREPARATION:**

Before starting any work, the containers must be completely clean and free of traces of wax or grease, so that the chemical union between the resins may be perfect.

We recommend a room temperature for mixing and application between 15°C and 30°C.

**DOSING AND PLASTICIZATION:**

Recommended ratio: powder/liquid 2:1 (14g of powder for 7ml of liquid).

Place the necessary amount of liquid in the glass, silicon or Dappen jar, in accordance with the amount of acrylic. Over the liquid, place the resin until the two are mixed together without having excess powder or liquid, if there is, dispose it and homogenize carefully with a spatula, avoiding the incorporation of air bubbles. Observe so that the mixture is liquid, however consistent.

Do not use the resin, when it stops sliding on the spatula (5 - 6 minutes after dosing, at 23°C).

Avoid touching with the fingers on the repaired surface.

**POLYMERIZATION:**

For polymerization, take the resin to a polymerizer, place the necessary amount of water, apply 20 pounds of compressed air and allow it to complete the polymerization, for 20 minutes.

Wait for total polymerization to start finishing and polishing.

For repairs, additions and rebases, proceed in the conventional way.

This resin can be used for orthodontic works.

High stability of colors and abrasion resistance.

NOTE - VIPI recommends for this product, solely mechanical polishing.

**FINISHING AND POLISHING:**

Finishing is started with bits and mills, passes to abrasive stones and is concluded with fine textured sandpapers and finishing Rods.

Polishing is performed with soft brush and polishing Rods.

The use of monomer with Crosslink will make polishing and brilliance easier, providing greater durability of the denture.

**Characteristics/Benefits:**

Rigorous selection and raw material receipt tests combined with successive tests during all manufacturing cycles, under ISO standard requirements, assure the specified results in the final product.

**Product Configurations:**

POWDER: Bottles of 20g, 80g, 225g, 450g, 1kg, 2,250kg and 5kg. Barrel of 25Kg and 50Kg.

LIQUID: Bottles of 50ml, 120ml, 250ml, 500ml and 1000 ml. 5L in can Gallon.

**Shades:**

Colorless, Pink and Black (only 80g bottle)

## Acrylic Resin

### Composition:

POWDER		LIQUID
Polymethylmethacrylate	Methyl methacrylate monomer	
Polypropylene		DMT
Pigments		Inhibitor

### Shelf Life:

POWDER: 10 years as of the date of manufacture.  
LIQUID: 2 years as of the date of manufacture.

### Infection control procedure

#### INSTRUCTIONS FOR DENTURE DISINFECTION

Note: *Instructions are not intended for routine denture care.*

	Method	Warning
<b>Cleaning &amp; Disinfection: Manual</b>	<ul style="list-style-type: none"><li>• Saturate a lint-free wipe with an alcohol quaternary-based solution.</li><li>• Wipe the entire appliance ensuring to wipe all grooves and edges well.</li><li>• Spray the entire surface with the alcohol quaternary-based solution to thoroughly wet the appliance - <i>maintain wet for the contact time as recommended per the manufacturer's instructions.</i></li><li>• After the manufacturer recommended contact time has elapsed, rinse the appliance under potable water for a minimum of 30 seconds.</li><li>• Allow the appliance to air dry completely.</li></ul>	<b>*Follow the manufacturer's instructions for solution contact time.</b>

### Recommendations to patients:

To increase the durability of the prosthesis, daily hygiene is recommended without the use of aggressive substances, such as acids or alkalis.

### ANVISA Registration n°:

POWDER: 10216040031  
LIQUID: 10216040032

### Technical Responsible:

Desiree M. Cortez - CRF-SP-41566

Check the correlation of the version of the instruction for use with the purchased product indicated on the packaging. To obtain the Instruction for Use in printed form at no additional cost, ask the customer service (SAC) by calling 0800 771 2226 (Brazil only) or sending an email to [atendimento@dentsplysirona.com](mailto:atendimento@dentsplysirona.com)

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