Resina Acrílica

Indicação:

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Acrylic Resin

Indication:

Elaboration of individual molding trays and test bases.

Characteristics/Benefits:

- Self-polymerizable resin of excellent dimensional stability.
- It allows excellent cost and benefit, considering the precision in molding.
- Rigorous selection and raw material receipt tests combined with successive tests during the entire manufacturing cycle, under ISO standard requirements, assure the specified results in the final product.
- Its composition allows a greater working time.
- Ease of handling and application in execution of denture work.

Presentations:

Bottles of 225 g and 2250 g.

Colors:

Ivory and Blue.

Composition:

POWDER

Polymethylmethacrylate

Benzoyl Peroxide

Biocompatible Pigments

Instructions for Use:

Shake the bottle before using.

PREPARATION:

Before starting any work, the parts to be used should be completely clean and free from traces of wax or fat, so that the chemical union between the resins may be perfect.

We recommend a room temperature for mixing and application between 15°C and 30°C..

DOSING AND PLASTICIZATION:

Recommended ratio: powder/liquid 2:1 (14g of powder for 7ml of liquid).

Place the necessary amount of fluid (VIPI FLASH) in the glass, silicon or Dappen jar, in accordance with the amount of acrylic. Over the liquid, place the resin until the two are mixed together without having excess powder or fluid, if there is, dispose it and homogenize carefully with a spatula, avoiding the incorporation of air bubbles. Observe so that the mixture is fluid, however consistent.

Do not use the resin, when it stops sliding on the spatula (5 - 6 minutes after dosing, at 23°C).

Avoid touching with the fingers on the repaired surface.

POLYMERIZATION:

For polymerization, take the resin to a polymerizer, place the necessary amount of water, apply 20 pounds of compressed air and allow it to complete the polymerization, for 20 minutes.

Wait for total polymerization to start finishing and polishing.

For repairs, additions and rebases, proceed in the conventional way.

This resin can be used for orthodontic works, however, we recommend the use of ORTO COR resin, especially formulated for such purpose.

High stability of colors and abrasion resistance.

NOTE - VIPI recommends for this product, solely mechanical polishing.

FINISHING AND POLISHING:

Finishing is started with bits and mills, passes to abrasive stones and is concluded with fine textured sandpapers and VIPI BRIL ACABAMENTO.

Polishing is performed with soft brush and VIPI BRIL POLIMENTO.

Warnings / Precautions:

- Follow the usual procedures for microbiological control, assembly, articulation and maintenance of the denture, as well as for individual protection.
- Exclusively for dental use.
- It can form flammable air/powder mixtures.
- As in all acrylic resins, handling of the product is recommended in ventilated places, using nitrile gloves and protective goggles; however, handling is not recommendable for pregnant women, nursing mothers and people with known allergies to this type of product.
- Conserve by protecting from the light and heat.
- Do not dispose the wastes in sewers and watercourses or in household garbage.

Contraindications:

There are no known adverse side effects, when the product is used according to usual procedures.

Acrylic Resin

Shelf Life:

10 years as of the date of manufacture.

ANVISA Registration nº.:

10216040031

Technical Responsible:

Luiz Paulo Ribeiro Junior - CRQ IV (Regional Council of Chemistry) 03251745

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