Silicone

Indicação:

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Silicon

Indication:

Setting silicon by condensation, of low viscosity and high hardness, especially elaborated for manufacture of denture walls.

Characteristics/Benefits:

- For use in conventional and microwave techniques.
- It involves innovative technique and easy application.
- It propitiates high fidelity in details of the model.
- It exempts the denture from plaster incrustations.
- It allows uniform and controlled polymerization.
- Economical and high yielding product.
- Best cost /benefit in the market, considering saving time in denture finishing.

Presentations:

KIT WITH:

- 1 flexible tube of 40 ml ACTIVATOR
- 1 flexible tube of 300 ml VIPI SIL BASE
- 1 manipulation jar
- 2 flexible spatulas for application and instructions manual

Observation: The ACTIVATOR and the VIPI SIL BASE are also sold separately.

Composition:

ACTIVATOR

Inorganic Components
Silicon Dioxide
Polydimethylsiloxane
Hydrocarbons
Red Food Coloring

VIPI SIL BASE

Dimethyl Polysiloxane
Silicon Dioxide
Biocompatible Pigments

Instructions for Use:

PREPARATION: The wax denture model is prepared conventionally and included in the muffle, where the plaster parts must be insulated normally. Retentions and characterizations in the teeth must be performed before waxing because the thickness of the VIPI SIL does not allow the removal and replacement of the teeth on the wall.

DOSING: To elaborate a normal size denture, 15 g of the BASE and 15 cm of the ACTIVATOR catalyst are necessary. The 15g of BASE are measured at the first mark of the measuring jar. Performing 1 round of the ACTIVATOR catalyzing paste result in 15 cm.

Considering the room temperature at 22°C, at this proportion we will have between 7 to 10 minutes. This time can be anticipated by increasing the amount of ACTIVATOR catalyst. Any alteration (between limits) in the dosing does not affect the properties of the final product. The only effect is a shorter (more ACTIVATOR catalyst) or more extended (less ACTIVATOR catalyst) pressing time.

APPLICATION: The VIPI SIL silicon must be mixed and applied with the spatula initially on the interdental gaps, in thin layers to avoid incorporation of air bubbles and later in the other areas of the denture, in a layer with thickness of 1 to 2 mm.

At the end of the work, the consistency of the mass increases allowing application of points of retention. We call the attention to the basic method in order to create retentions.

APPLICATION OF RETENTIONS: Have at hand, before starting the work, pieces of silicon used in previous works, with the approximate size of 1.5 cm, in sufficient number to cover the model. After the total coverage of the model with VIPI SIL, apply the pieces of silicon on the surface, assuring an efficient retention.

INCLUSION: After final curing close the muffle and complete with plaster.

The VIPI SIL silicon can be used both in the microwave cooking technique and in the conventional technique.

Observations: The presence of Silicon fluid when opening the Flexible Tube for the first time is normal. When using the Flexible Tube of the ACTIVATOR catalyst, close it immediately to avoid the entry of air moisture.

Warnings/Precautions:

- Follow the usual procedures for microbiological control, assembly, articulation and maintenance of the denture, as well as for individual protection.
- Exclusively for dental use;
- Flammable product;
- Maintain between 5°C and 27°C;
- It can cause irritability;
- This resin's handling is not recommendable to pregnant women, nursing mothers and people with known sensitivities to this type of product.
- Keep away from light and hear;
- Do not dispose the wastes in sewers and watercourses or in household garbage;

Silicon

Shelf Life:

2 years of the date of manufacture.

Product exempted from registration by the Ministry of Health.

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