### Opacificador

#### Indicação:

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### Opacifier

#### Indication:

Metal reinforcements of total denture, metal saddles of removable partial denture, protocol bars over implants and on any metal surface to be coated with acrylic resin.

#### Contra indications:

There are no known adverse side effects, when the product is used according to usual procedures.

#### Warnings/Precautions:

- Follow the usual procedures for microbiological control, assembly, articulation and maintenance of the denture, as well as for individual protection.
- Exclusively for dental use.
- FLUID: Flammable, non-toxic material (Flash point: 10°C/50°F).
- RESIN: It can form flammable air/powder mixtures.
- As in all acrylic resins, handling of the product is recommended in ventilated places, using nitrile gloves and protective goggles; however, handling is not recommendable for pregnant women, nursing mothers and people with known sensitivities to this type of product.
- Conserve by protecting from the light and heat.
- Do not dispose the wastes in sewers and watercourses or in household garbage.

#### Instructions for Use:

After adequate machining of the metal structure, with the proper mechanical retentions, apply a sandblast and then wash and dry well. To improve the adherence of the opaque over the metal, application of a thin layer of VIPI TONE SILANO is recommended. After drying for 5 minutes in open air, expose the structure with SILANO for another 10 minutes to a source of dry heat, between 80°C and 100°C. According to the technique to be used, place a few drops of the fluid (Photopolymerizable STG Fluid or Self-polymerizable STG Fluid) in a Dappen jar. Shake the content of the bottle of the OPACIFIER powder, in the chosen color, and add a portion on that fluid, sufficient to form a well fluid mixture. Apply with dry brush, thin layers over the metal surface, until the entire metal is covered, as many coatings as needed.

If the fluid used in the mixture is the photopolymerizable, subject the structure with the applied opaque, to a photopolymerizer device, until the opaque is well hardened.

The characteristics of the device's polymerization must be taken into account. As reference we indicate: UNIX (Kulzer) 180s / Fotoceram 300s / Xenon pro 180s / HI-LIGHT180s / EDGLUX 240s

If the fluid used is self-polymerizable, after application, allow it to dry for 20 minutes.

#### Characteristics/Benefits:

- Ease of handling.
- Aesthetic and effective results in the opacification of any metal structure.
- Independence in the mode of polymerization (photopolymerization or self-polymerization).
- A union agent (SILANE) is available to prevent infiltrations between the resin and the metal.
- It does not need special solvents for dilution and cleaning.
- It has photopolymerizable and self-polymerizable fluid.

#### Presentations:

KIT WITH:

3 bottles of VIPI TONE OPACIFIER of 20 g each;

1 Photopolymerizable Fluid of 15 ml;

1 Self-polymerizable Fluid of 15 ml;

1 Union Agent VIPI TONE SILANO 5 ml

Observation: All items of the KIT are also sold separately for replacement.

#### Colors:

VIPI TONE OPACIFIER: Light Pink, Dark Pink, Dentin 62

#### Composition:

OPACIFIER (POWDER)

Polymethylmethacrylate

Titanium Dioxide

Biocompatible Pigments

SILANE

Alcohol

Methacrylsilane

SELF POLYMERIZABLE STG FLUID

Methylmethacrylate

DMT

EDMA (Crosslink)

Inhibitor

PHOTOPOLYMERIZABLE STG

FLUID

BDMA

Initiator

Inhibitor

#### Shelf Life:

POWDER: 10 years as of the date of manufacture.

FLUID: 2 years as of the date of manufacture.

## Opacifier

**ANVISA Registration n°:**

SILANE AND OPACIFIER: 10216040014

SELF-POLYMERIZABLE AND PHOTOPOLYMERIZABLE FLUIDS: 10216040032

**Technical Responsible:**

Luiz Paulo Ribeiro Junior – CRQ IV (Regional Council of Chemistry) 03251745

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