### Resina

#### Indicação:

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### Resin

#### Indication:

Red photopolymerizable resin, specially prepared for the elaboration of caps, elaboration of permanent dentures and dentures on implants; attachment contour of telescopic precision and conic crowns; Contour of pillars of individual implants; Contour of adhesive bridges, inlays and onlays; Contour of crowns and bridges for pantographic transmission; Intraradicular pin [filling and capture structure]; Customized molding of implant protractors; Locking of multiple drag molding on implants; Coping for elaboration of permanent dentures; Caps for permanent denture; Structures for pressed ceramic like inlay laminates; direct casting RPD (Removable Partial Denture), without duplicate model.

#### Characteristics / Benefits:

- Direct application on the model
- Viscous consistency [thixotropic], like gel
- Ideal consistency, immediate use
- Easy and fast to mold
- Photopolymerize in thin or continuous layers
- Flexibility, fast preparation for casting models
- Fast curing
- The surface can be smoothed with electric probe
- It does not leave residues to the ring after burning

#### PHYSICAL PROPERTIES:

- Bending resistance: Approximately 50 MPa
- Elasticity module: Approximately 1200 MPa
- Appearance, consistency, viscosity: Bubble-free, firm gel 150 – 350Pas
- Vickers Hardness: over 12 HV
- Color: Red
- Curing Depth: > 0.5 cm
- Residue-free, ash-free burning: Yes

#### Presentations:

1 syringe with 3gr and 2 measuring nozzle

#### Color:

Red

#### Composition:

Methacrylates  
Activators

#### Instructions for Use:

Apply the product with the metal nozzle over the area where the cap will be built. Take to the photopolymerizing oven and perform the application in thin layers, avoiding the non-photopolymerization of the material. The whole application can be done by layers with photopolymerizing device for 20 seconds.

In case you perform the resin application with a brush, place a small portion on a glass plate and, with small portions apply it spreading on the entire prepared area.

**PATTERN DIFFERENCE:** In the powder / fluid brushing technique the process is slow, due to the need of waiting for resin polymerization, without having in mind the influence of proportion.

#### Warnings / Precautions:

- Exclusively for dental use.
- Keep the syringe closed after use.
- To prevent resin photopolymerizing in to the metal nozzle, keep product protected from the light.
- Protect from the light and heat.
- Store between 2°C and 28°C.
- Do not dispose the wastes in sewers and watercourses or in household garbage.

**Shelf Life:** 4 years from manufacturing date.

**ANVISA Registration n°:** 10216040024

**Technical Responsible:** Hélio Wiebeck - CRQ (Regional Council of Chemistry) 04302790

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**Madespa S.A.**