Document Detail

Type: 1100-PIR-DFU
Document No.: 1100-PIR-DFU 0035[7]
Title: FICHA TÉCNICA VIPIBLOCK PMMA
Comment:
Status: CURRENT
Effective Date: 24-Jul-2023

Approval

<u>Owner Role</u>	<u>Sign-off By</u>	<u>Sign-off Date</u>	
1100-Research and Development Brazil Research and Development	Luiz Crepaldi	24-Jul-2023 6:34 pm	GMT
1100-Marketing Approver Marketing Approver	Fabio Domingues	24-Jul-2023 5:59 pm	GMT
1100-Regulatory Approver Brazil Regulatory Approver	Ana Carolina Nilson	24-Jul-2023 5:59 pm	GMT
1100-Document Control Document Control	Priscilla Silva	24-Jul-2023 5:58 pm	GMT
1100-PIR-ENGENHARIA Aprovador Brasil Pirassununga Engenharia Aprovador	Bruno Conti	24-Jul-2023 5:54 pm	GMT
1100-PIR-QAAR Aprovador Brasil Pirassununga QAAR Aprovador	Gabriel Ocanha	24-Jul-2023 5:55 pm	GMT
1100-Clinical Education Appro Clinical Education Approver	Ana Zimbaldi	24-Jul-2023 5:56 pm	GMT

Reference

<u>Document No.</u>	<u>Title</u>	<u>Relation</u>
1100-PIR-DMR 0001 [19]	Multilayer PMMA Disc	Related

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ENGLISH

Indication:

VIPI BLOCK TRILUX and VIPI BLOCK MONOCOLOR:
PMMA block for temporary (single or multiple) crowns and bridges, created and produced by CAD/CAM system.

Contraindications

- Contraindicated for people allergic to any of the components

Warnings and precautions:

- Use only by dental professionals trained in operating CAD/CAM systems.
- Avoid overheating during machining and finishing.
- Dust created during machining can cause irritation to the eyes, skin and respiratory tract. Use of personal protective equipment (mask, glasses) is recommended.
- The structure must be cleaned and disinfected before being placed in the patient's mouth.

Recommendations to the patient:

- To increase the durability of the prosthesis, daily cleaning without the use of aggressive substances, such as acids and alkalis, is recommended.
- In case of accidental ingestion or allergic reactions to the medical device it is recommended to visit a dentist or doctor as soon as possible.

Disposal: The product does not require special conditions for disposal, however, it is recommended to consult the legislation in force in the country before proceeding with final disposal.

Additional information:

- Check the correlation of the version of the instruction for use with the purchased product indicated on the packaging.
- To obtain the Instruction for Use in printed form at no additional cost, ask the customer service (SAC) by calling 0800 771 2226 (Brazil only) or sending an email to atendimento@dentsplysirona.com.
- Do not use if the packaging is damaged.
- FORBIDDEN REPROCESS

Instruction for Use:

VIPI BLOCK is used in the machining technique to make crowns and bridges, such as long-lasting, cemented or screwed temporary crowns and bridges (Trilux and Monocolor) machined by CAM System from different manufacturers.

1. Attach the block to the machine following the equipment manufacturer's instructions.
2. When using Trilux Block, the color variation in the element is given by the positioning of the piece according to the horizontal layers of colors of the block, according to figures 1, 2 and 3 below.
Check the availability in changing the vertical position in the CAM software used.
3. After machining, remove the block from the equipment and remove the element from the block.
4. To finish the finished work, use cutters, sandpaper and rubbers.


Limits Posterior:

Trilux and Monocolor Blocks: Maximum 2 fixed elements plus one suspended element.

Minimum thickness: 1 mm occlusal and 0.6 mm cervical.

Minimum connection diameters: Anterior: 3.6 mm; Rear: 4.4 mm.

The temporary block can be fixed with all cements for provisional fixation. Temporary cements with Eugenol can decrease the effectiveness of resin cements.

Images showing the element positioned at different heights of the block, changing the color of the incisal edge.
For each CAD/CAM system there is a different form of positioning.

Figure 1: Piece positioned in the upper region of the block: More incisal.

Figure 2: Piece positioned in the lower region of the block: Less incisal.

Figure 3: By tilting the part, it is possible to smooth the gradient.

Aesthetic suggestion table:

COLOR TYPES	COLORS	INDICATION SUGGESTED
TRILUX	A1, A2, A3, A3,5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4, HW2	For crowns and bridges (single or multiple).
MONOCOLOR	A1, A1T, A2, A3, A3,5, A4, B1, B1T, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4, HW2	For crowns and bridges (single or multiple).

Performance characteristics:

The product meets the applicable requirements of ISO 10477: 2020 Dentistry - Polymer-based crown and veneering materials.

Symbols on VIPI BLOCK TRILUX and VIPI BLOCK MONOCOLOR labels:

Manufacturer	Expiration Date
Medical Device	Manufacture Date
Reference Number	Consult Instructions For Use
Batch or Lot Number	Keep Dry
Caution	Keep Away From Sunlight

Composition:

Methyl Polymethacrylate, Pigments, EDMA and Fluorescent.

Expiration date:

10 years after the manufacturing date.

Packaging:

6 pieces of VBS model

1 piece on other models

Storage:

keep in a dry place and protected from sunlight.

Technical Responsible: Desiree M. Cortez - CRF-SP-41566 (information applicable to Brazil)

VIPI BLOCK TRILUX and VIPI BLOCK MONOCOLOR - ANVISA No. 10216040029 (Resinas Acrílicas Odontológicas)

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